

Marco Rocchini INFN - Istituto Nazionale di Fisica Nucleare FIRENZE DIVISION

SPIDER: the Silicon Ple DEtectoR



Outline

- Low-energy Coulomb excitation
- Shape of ground and excited states via Coulomb excitation
- What is **SPIDER**?
 - Energy resolution, modularity, simulations
- Pulse height correlations
- Radiation damage
- Online RBS analysis with SPIDER
- The GALILEO campaign
- The AGATA campaign



measurements



Low-Energy Coulomb Excitation

SPIDER

Pulse Height Correlations

Radiation Damage


Online Target Analysis


The GALILEO Campaign

The AGATA Campaign

Inelastic scattering between two interacting nuclei, in the "purely" electromagnetic regime









Low-Energy Coulomb Excitation

SPIDER


Pulse Height Correlations


Radiation Damage

Online Target Analysis

The GALILEO Campaign

- Inelastic scattering between two interacting nuclei, in the "purely" electromagnetic regime
- The time-dependent electromagnetic field between the two nuclei can induce excitations







Low-Energy Coulomb Excitation

SPIDER

Pulse Height Correlations

Radiation Damage


Online Target Analysis


The GALILEO Campaign


The AGATA Campaign

- Inelastic scattering between two interacting nuclei, in the "purely" electromagnetic regime
- The time-dependent electromagnetic field between the two nuclei can induce excitations
- The nuclei then de-excite; in Coulex we are mostly interested in γ-ray emission

 $d\sigma_{clx}$







 $\frac{d\sigma_{Ruth}}{} \cdot P\left(i \longrightarrow f\right)$

$$P\left(0_{1}^{+} \longrightarrow 2_{1}^{+}\right) = F\left(\theta, E_{P}\right) B(E2) \left[1 + 1.32 \frac{A_{P}}{Z_{T}} \frac{\Delta E}{\left(1 + \frac{A_{P}}{A_{T}}\right)} Q_{s}\left(2^{+}\right) K\left(\theta, E_{P}\right)\right]$$



Low-Energy Coulomb Excitation

SPIDER

Pulse Height Correlations


Radiation Damage


Online Target Analysis


The GALILEO Campaign

The AGATA Campaign

- Inelastic scattering between two interacting nuclei, in the "purely" electromagnetic regime
- The time-dependent electromagnetic field between the two nuclei can induce excitations
- The nuclei then de-excite; in Coulex we are mostly interested in γ-ray emission







$$\frac{d\sigma_{clx}}{d\Omega} = \frac{d\sigma_{Ruth}}{d\Omega} \cdot P(i \longrightarrow f)$$

Example: first 2+ state in an even-even target nucleus

$$P\left(0_{1}^{+} \longrightarrow 2_{1}^{+}\right) = F\left(\theta, E_{P}\right) B(E2) \left[1 + 1.32 \frac{A_{P}}{Z_{T}} \frac{\Delta E}{\left(1 + \frac{A_{P}}{A_{T}}\right)} Q_{s}\left(2^{+}\right) K\left(\theta, E_{P}\right)\right]$$

Access to: transition probabilities, spectroscopic quadrupole moments



Quadrupole Sum Rules

Low-Energy Coulomb Excitation

SPIDER

Pulse Height Correlations


Radiation Damage

Online Target Analysis

The GALILEO Campaign

The AGATA Campaign

- ▶ Unique feature of low-energy Coulomb excitation ⇒ Possible to get relative signs of transitional matrix elements, and spectroscopic quadrupole moments of excited states with their sign
- Quadrupole Sum Rules \Rightarrow (β, γ) deformation parameters for g.s. and excited states in a model-independent way



K. Kumar, Phys. Rev. Lett. 28 (1972) 249 & D. Cline, Annu. Rev. Nucl. Part. Sci. 36 (1986) 683



Quadrupole Sum Rules

Low-Energy Coulomb Excitation

SPIDER

Pulse Height Correlations

Radiation Damage

Online Target Analysis


The GALILEO Campaign


The AGATA Campaign

- ► Unique feature of low-energy Coulomb excitation ⇒ Possible to get relative signs of transitional matrix elements, and spectroscopic quadrupole moments of excited states with their sign
- Quadrupole Sum Rules \Rightarrow (β, γ) deformation parameters for g.s. and excited states in a model-independent way


$$\langle i | Q^{2} | i \rangle = \frac{\sqrt{5}}{\sqrt{2I_{i}+1}} \sum_{t} \langle i | | E2 | | t \rangle \langle t | | E2 | | i \rangle \begin{cases} 2 & 2 & 0 \\ I_{i} & I_{i} & I_{t} \end{cases}$$

$$\langle i | Q^{3} \cos(3\delta) | i \rangle = -\frac{\sqrt{35}}{\sqrt{2}} \frac{1}{2I_{i}+1} \sum_{tu} \langle i | | E2 | | t \rangle \langle t | | E2 | | u \rangle \langle u | | E2 | | i \rangle \begin{cases} 2 & 2 & 2 \\ I_{i} & I_{t} & I_{u} \end{cases}$$









equivalent to

 (β, γ)

K. Kumar, Phys. Rev. Lett. 28 (1972) 249 & D. Cline, Annu. Rev. Nucl. Part. Sci. 36 (1986) 683



Excited 0+ States in Mid-Mass Nuclei

Low-Energy Coulomb Excitation

SPIDER

Pulse Height Correlations

Radiation Damage

Online Target Analysis

The GALILEO Campaign





Excited 0+ States in Mid-Mass Nuclei

Low-Energy Coulomb Excitation


SPIDER

Pulse Height Correlations

Radiation Damage

Online Target Analysis

The GALILEO Campaign





Experimental Setup

Low-Energy Coulomb Excitation

SPIDER


Pulse Height Correlations

Radiation Damage

Online Target Analysis

The GALILEO Campaign

- γ-ray detector:
 - To measure γ-ray yields for the different transitions
- Particle detector:
 - To select the events of interest
 - To distinguish between projectile and target nuclei
 - To perform the Doppler correction
 - To select the scattering angle
 - To monitor the experiment





SPIDER: the Silicon Ple DEtectoR

Low-Energy Coulomb Excitation

SPIDER


Pulse Height Correlations

Radiation Damage

Online Target Analysis

The GALILEO Campaign

- Modular segmented silicon detector, originally designed for low-energy Coulomb-excitation measurements
- Independent sectors, 8 strips + guard ring
- Detector thickness ~ 300 μm, dead layers ~ 50 nm in the junction (front) side and ~ 350 nm in the ohmic (rear) side
- Cone configuration (7 sectors) at backward angles: 8.5 cm from the target $\Rightarrow \Delta\Theta = 37.4^{\circ}$, $\Omega/4\pi = 17.3\%$









Energy Resolution

Low-Energy Coulomb Excitation

SPIDER

Pulse Height Correlations


Radiation Damage


Online Target Analysis

The GALILEO Campaign

The AGATA Campaign

• nTD silicon detectors, high quality, high energy resolution (less than ~ 1% for 5.5-MeV α particles)





M. Rocchini, K. Hadyńska-Klęk, A. Nannini et al., NIMA 971 (2020) 164030



Modularity

Low-Energy Coulomb Excitation

SPIDER

Pulse Height Correlations


Radiation Damage

Online Target Analysis

The GALILEO Campaign

The AGATA Campaign


The detector can be mounted in different arrangements, using from 1 to 8 sectors, giving different Doppler-correction capabilities and angular coverages



Polar and solid angles with an 8.5-cm distance from the target:

Strip	$ heta_{min}$	$ heta_{mean}$	$ heta_{max}$	${\it \Delta} heta$	arOmega
	[deg]	[deg]	[deg]	[deg]	[srad]
7	123.5	125.4	127.5	4.0	0.046
6	127.5	129.6	131.8	4.3	0.046
5	131.8	134.0	136.4	4.6	0.045
4	136.4	138.7	141.2	4.8	0.043
3	141.2	143.6	146.1	5.0	0.040
2	146.1	148.6	151.2	5.1	0.035
1	151.2	153.7	156.3	5.1	0.030
0	156.3	158.8	161.3	5.1	0.028

Simulated ⁶⁰Ni on ²⁰⁸Pb at the safe energy



M. Rocchini, K. Hadyńska-Klęk, A. Nannini et al., NIMA 971 (2020) 164030



Geant4 Simulations

Low-Energy Coulomb Excitation

SPIDER


Pulse Height Correlations

Radiation Damage


Online Target Analysis


The GALILEO Campaign

- Geant4 package available to simulate in detail each experiment
- Input from GOSIA (reference code for Coulomb-excitation cross-section calculations)
- Possibility to include both the beam and the target excitations



M. Rocchini, K. Hadyńska-Klęk, A. Nannini et al., NIMA 971 (2020) 164030







Pulse Height Correlations

Low-Energy Coulomb Excitation

SPIDER

Pulse Height Correlations


Radiation Damage

Online Target Analysis

The GALILEO Campaign

The AGATA Campaign

- Time coincidences between SPIDER signals (rejected in the analysis)
 - A = double hits (~ 0.2% of the total statistics)
 - B,C = real signals in coinc. with noise and lowamplitude signals
 - D = noise coincidences



M. Rocchini, K. Hadyńska-Klęk, A. Nannini et al., NIMA 971 (2020) 164030



Pulse Height Correlations

Low-Energy Coulomb Excitation

SPIDER


Pulse Height Correlations

Radiation Damage


Online Target Analysis

The GALILEO Campaign

- Time coincidences between SPIDER signals (rejected in the analysis)
 - A = double hits (~ 0.2% of the total statistics)
 - B,C = real signals in coinc. with noise and lowamplitude signals
 - D = noise coincidences
- Negative polarity signals









Radiation Damage

Low-Energy Coulomb Excitation

SPIDER

Pulse Height Correlations


Radiation Damage

Online Target Analysis

The GALILEO Campaign

The AGATA Campaign

Permanent lattice defects or temporary effects (temperature, auto-recovered displacements ...)



M. Rocchini, K. Hadyńska-Klęk, A. Nannini et al., NIMA 971 (2020) 164030



Online Target Analysis

The high energy resolution of SPIDER allows for the online application of the Rutherford-backscattering technique (e.g., using the SIMNRA code)

- Target thickness, contaminants, carbon-built-up layer, etc...
- Example in 66Zn on 208Pb with GALILEO:
 - > 208Pb thickness = 616(11) mg/cm², ¹²C built-up = 114(12) mg/cm²



Pulse Height

Correlations

Low-Energy

Coulomb


Excitation


SPIDER

Damage

Online Target Analysis

The GALILEO Campaign








Online Target Analysis

The high energy resolution of SPIDER allows for the online application of the Rutherford-backscattering technique (e.g., using the SIMNRA code)

Target thickness, contaminants, carbon-built-up layer, etc...



Pulse Height

Correlations

Low-Energy

Coulomb


Excitation

SPIDER

Radiation Damage

Online Target Analysis

The GALILEO Campaign





The Silicon Ple DEtectoR

GALILEO

Low-Energy Coulomb Excitation

SPIDER

Pulse Height Correlations


Radiation Damage

Online Target Analysis

The GALILEO Campaign

- 25 HPGe Compton-suppressed detectors (GASP type)
- FWHM (@1332.5 keV) < 2.4 keV
- Efficiency (@1332.5 keV) = 2.1%
- Full digital electronics (takes advantage of the developments made for AGATA):
 - Trigger-less mode
 - Typical operational rate ~ 20 kHz/det
 - Common clock synchronization
 - Local data processing





GALILEO with SPIDER (2016 - 2019)





INFN

SPIDER


Pulse Height Correlations

Radiation Damage

Online Target Analysis

The GALILEO Campaign

The AGATA Campaign



1000 1200 1400

014311

800

Energy [keV]

200

400



Type-II Shell Evolution in 94Zr

Low-Energy Coulomb Excitation


SPIDER

Pulse Height Correlations

Radiation Damage

Online Target Analysis

The GALILEO Campaign





Shapes in ¹³⁰Xe for ¹³⁰Te 0vββ Decay

Low-Energy Coulomb Excitation


SPIDER

Pulse Height Correlations

Radiation Damage

Online Target Analysis

The GALILEO Campaign





Coulex of 58Ni and 116Sn

Low-Energy Coulomb Excitation


SPIDER

Pulse Height Correlations

Radiation Damage

Online Target Analysis

The GALILEO Campaign





The Silicon Ple DEtectoR

AGATA



Low-Energy Coulomb Excitation

SPIDER


Pulse Height Correlations


Radiation Damage


Online Target Analysis

The GALILEO Campaign

- Last generation γ-ray spectrometer composed of high-segmented HPGe detectors
- Employs advanced PSA and γ-ray tracking methods to avoid Compton-suppressors and garantee high efficiency







AGATA with SPIDER (2022 - ongoing)





AGATA with SPIDER (2022 - ongoing)




SPIDER

Pulse Height Correlations

Radiation Damage

Online Target Analysis

The GALILEO Campaign





AGATA+SPIDER 2nd and 3rd Campaign

Low-Energy Coulomb Excitation


SPIDER

Pulse Height Correlations

Radiation Damage

Online Target Analysis

The GALILEO Campaign






Summary

SPIDER is a segment, modular silicon detector designed for low-energy Coulomb excitation measurements with modern γ-ray spectrometers

The detector is characterized by an excellent energy resolution, which makes it possible to take into account the target thickness in the Doppler correction of the γ-ray energy

We studied radiation damage and pulse height correlations for our applications

The detector is "as simple as possible", making it a very reliable instrument (100% success rate)



measurements



Marco Rocchini INFN - Istituto Nazionale di Fisica Nucleare FIRENZE DIVISION

THANK YOU FOR THE ATTENTION



